

# Specifications for organiQ THERAPY BMR

## BASIC CHARACTERISTICS

The device model	OrganiQ THERAPY BMR
Number of programs	11+4
Duration of therapy	15 minutes
Frequency magnetic field	1,2 - 9999 Hz
Strength of the magnetic field	0.012 - 50 Gauss
The voltage of the input power	110 - 240 VAC
The frequency of the input power	47 - 63 Hz
Power input	42 VA
Device dimensions (approx.)	8.9x2.8x7.8 inches
GVWR (approx.)	1.98 lbs

## OPERATING CONDITIONS AT WORK

Temperature	5°F to 104°F
Relative humidity	30% to 75%
Atmospheric pressure	20.6 inHg to 31.3 inHg

## Reference centers



Available soon in Canada, Qatar, Pakistan, India, Australia and Japan.

## Contact us

[www.hsttsolutions.com](http://www.hsttsolutions.com)



HSTT Solutions  
8588 Katy Freeway Suite 230  
Houston, Texas 77024

100%  
NONINVASIVE

organiQ  
THERAPY

# Bio-Magnetic Resonance



# What is organiQ THERAPY BMR?

**organiQ THERAPY** Bio-Magnetic Resonance (BMR) is an innovative device that provides pulsating bio-magnetic resonance therapy. It was developed based on research which respects the physiological processes in the body.


**organiQ THERAPY BMR** excites the movement of electrically charged particles (ions) in the cell using a pulsed magnetic field and a specific frequency, so that there is several times higher permeability of the cell membrane. This improves the input of nutrients and oxygen to cells, and the output of toxins and heavy metals.

**organiQ THERAPY BMR** offers a choice of 11 therapies and 4 programs to stimulate brain waves (alpha, beta, delta and theta). Duration of one session is 15 minutes.

**organiQ THERAPY BMR** provides a completely natural, painless and non-invasive form of therapy intended for all age groups. Differences compared to the devices used in the field of magnetic resonance:

- ▶ **organiQ THERAPY BMR**, by selecting the program sets the frequency, the magnetic field density and signal formats. This means that these allowing functions can be targeted (according to doctor's orders) effectively to treat the specific type of pathology.
- ▶ **organiQ THERAPY BMR** energy is released in the form of light energy (photons) and does not cause micro damage of the soft tissue.

## Application



### PAIN CLINIC

Eliminating or significantly reducing the pain of all etiologies



### RHEUMATOLOGY

Rheumatoid arthritis, lumbosacral syndrome, periartrosis of shoulder joint, tennis elbow, sudek syndrome



### NEUROLOGY

Resolving insomnia and sleep problems



### ORTHOPEDY

Bone fractures, sports injuries, chronic pain (bone spurs, ...)



### PHYSICAL MEDICINE

Postoperative rehabilitation trauma



### SPORTS MEDICINE

All forms of mechanical injuries (muscles, tendons, bones, joints)

## Advantages

- ▶ A completely natural way of therapy
- ▶ Respects the physiological processes in the body
- ▶ Painless therapy, not invasive or aggressive
- ▶ Therapy for all age groups
- ▶ The treatment is performed in a lying or sitting position
- ▶ One treatment lasts 15 minutes
- ▶ One therapy is sufficient daily or every other day

## Contraindications

- ▶ Pacemaker
- ▶ Hyperactivity of the thyroid
- ▶ Epilepsy
- ▶ Extensive internal bleeding
- ▶ Shrapnel in the body
- ▶ Expectant mothers